Hydrogum - Orthoprint Tropicalgin - Neocolloid

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Alginate

Alginate dental impression material

Hydrogum

Orthoprint

Tropicalgin

Neocolloid

EN. Use protective goggles, filtering face mask suitable for powders (FFP2 or FFP3), gloves and work clothes.

PREPARING THE DEVICES: STEP-BY-STEP INSTRUCTIONS

A. 453g/500g BAG (Hydrogum/ Orthoprint/ Tropicalgin/ Tropicalgin Normal/ Tropicalgin Extra Fast/ Neocolloid)

1.1

1.2

1.3


1.4

1.5


1.6


1.7

1.8


1.9

1.10


B. 18g SINGLE DOSE BAG (Orthoprint/ Tropicalgin)

2.1


2.2


2.3

IMPRESSION TAKING PROCEDURE


3.1


3.2

DISINFECTING AND CLEANING THE DEVICE: GENERAL STEP-BY-STEP INSTRUCTIONS

4.1

4.2

EN – Hydrogum – Orthoprint – Tropicalgin – Tropicalgin Normal – Tropicalgin Extra Fast – Neocolloid

Instructions for use

1. INTENDED USE

Alginates dental impression material.

2. PRODUCT DESCRIPTION*

Hydrogum: an historical Zhermack alginate. Mint scent.

Orthoprint: alginate with short time in mouth, ideal for orthodontics. Vanilla scent.

Tropicalgin, Tropicalgin Normal and Tropicalgin Extra Fast: highly-versatile chromatic alginate. The chromatic variation accompanies the different processing phases: fuchsia during mixing, orange while working and positioning on impression tray, yellow during the time in mouth until curing. Lemon and tropical scent.

Neocolloid: the alginate with a longer time in mouth, recommended for removable prostheses. Mint and sage scent.

**The information above refers to the Zhermack alginate range.*

3. PACKAGING

453 g bag (Tropicalgin, Tropicalgin Normal, Tropicalgin Extra Fast)

500 g bag (Hydrogum, Orthoprint, Neocolloid)

18 g single-dose bag (Orthoprint, Tropicalgin)

4. COMPOSITION

Hydrogum: silica fillers, calcium sulphate, alginate, phosphate, potassium fluorotitanate, magnesium, hydrocarbons, zinc oxide, pigments, mint scent.

Orthoprint: silica fillers, calcium sulphate, alginate, phosphate, potassium fluorotitanate, magnesium, hydrocarbons, zinc oxide, pigments, vanilla scent.

Tropicalgin, Tropicalgin Normal and Tropicalgin Extra Fast: silica fillers, calcium sulphate, alginate, phosphate, potassium fluorotitanate, hydrocarbons, pH indicator, pigment, lemon and tropical scent.

Neocolloid: silica fillers, calcium sulphate, alginate, phosphate, potassium fluorotitanate, magnesium, hydrocarbons, zinc oxide, pigments, mint and sage scent.

5. COMPATIBLE DEVICES

- With Hydrogum, Orthoprint, Tropicalgin, Tropicalgin Normal, Tropicalgin Extra Fast and Neocolloid use the dosing set with a powder-to-water ratio of 9 g:18 ml.
- The material is compatible with study model or working model gypsums (type 3 gypsum e.g. Elite Model - Zhermack or type 4 gypsums e.g. Elite Rock - Zhermack).

6. INDICATIONS FOR USE

The devices are indicated to be used by dental practitioners (users authorized to operate on patients) in dental surgeries (public or private health facilities with the necessary legal authorizations) to take impressions on dental patients of any age.

7. CONTRAINDICATIONS

Do not use on patients with known hypersensitivity to any of the components. For further information, contact Zhermack.

8. SIDE EFFECTS

Irritation, redness or signs of hypersensitivity may occur in the case of allergy to any of the components. In case of hypersensitivity, discontinue use immediately and adopt measures to ensure patient safety.

9. CLINICAL BENEFITS

Impression material for accurate registrations of the dimensions of tissues and/or oral prosthetic components and the spatial relationships between them.

10. PRECAUTIONS

- There is insufficient data to support the use of alginate for impressions on exposed bone. The alginate is only intended for use on and around intact mucosa.
- Wear appropriate protective goggles, filtering face mask suitable for powders (FFP2 or FFP3), work clothes and gloves.
- All products must be used at room temperature. Allow the product to reach room temperature before use. High temperatures reduce (speed up) working and setting times, while low temperatures increase (slow down) these times.
- Do not use mixed product if excessively fluid.
- Do not use if the mixture is not homogeneous.
- In order to improve the retention between alginate and impression tray, it is advisable to use a perforated impression tray.
- Do not overfill the impression tray to reduce the gag reflex and risk of suffocation or ingestion.
- Cut out and keep the batch number and expiry date on the bag.
- Keep product accessories pack for the duration of use of the product.
- Do not use the material if the batch number and expiry date on the outer packaging are missing or illegible. Do not use after the expiry date.
- Check that the impression is intact and make sure there are no material residues to be removed from the patient's mouth.
- Use clean and properly disinfected and sterilised or ready-to-use accessories.
- Store in an airtight container used solely for alginate powder. Check that the container is clean and dry.
- Do not use the product if the packaging is damaged.
- The product has not been tested for compatibility with products other than those with which it is supplied. Using non-original accessories may have an adverse effect on the precision of the impression. Use of non-original accessories is therefore not recommended.
- Times for clinical use (see technical data table in paragraph 14) are understood to be from commencement of mixing with distilled water at 23°C/73°F. The temperature and hardness of the water may affect the times for clinical use of the material. Cold water retards curing; hot water accelerates it.
- Use distilled water. Changes in water hardness can affect setting times.

• Do not contaminate the material or the teeth/surfaces with the non-polymerized acrylic resins or bis-acrylic surface layer. Residues of these materials could interfere with product jellification. If the surfaces are contaminated, it is advised to clean them with ethyl alcohol or other solvents.

11. WARNINGS

- Do not use the material if the batch number and expiry date on the outer packaging are missing or illegible. Do not use after the expiry date.
- Avoid contact with the eyes. IN CASE OF CONTACT WITH EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. In the event of irritation, consult a doctor.
- The alginate contains respirable crystalline silica. Avoid prolonged or repeated inhalation of powders. Use a filtering face mask suitable for powders (FFP2 or FFP3). When using, do not eat, drink or smoke. In the event of inhalation and/or malaise, consult a doctor. Mix the water and powder with caution, if possible away from the patient.
- Do not swallow intentionally. To avoid unintentional ingestion caused by material jellification problems, comply with the information concerning use and storage provided in these instructions. In the event of ingestion and/or malaise, consult a doctor.
- To avoid damaging the periodontal tissue, check that the impression is intact and make sure there are no material residues to be removed from the patient's mouth. After removing the impression, have the patient rinse their mouth. Check there are no material residues in the gum area.

12. STEP-BY-STEP INSTRUCTIONS

12.1 PRELIMINARY OPERATIONS

Use protective goggles, filtering face mask suitable for powders (FFP2 or FFP3), gloves and work clothes.

Handle and dispense outside the surgical area in a clean place away from splashes, sprays or aerosols of body fluids.

In case of mechanical or automatic mixing, follow the instructions provided by the machine manufacturer.

12.2 MATERIAL PREPARATION

A. 453 g/500 g BAG (Hydrogum/ Orthoprint/ Tropicalgin/ Tropicalgin Normal/ Tropicalgin Extra Fast/ Neocolloid)

1. Open the bag (fig. 1.1) and pour the alginate into an airtight, clean and dry container (fig. 1.2).
2. Seal the container properly (fig. 1.3). Cut out the top part of the bag bearing the variable data (batch number, expiry date and UDI) (fig. 1.4) and affix it to the lid of the container (fig. 1.5).
3. Choose a suitable impression tray. Use the dosing set with a powder-to-water mixing ratio of 9 g:18 ml.
4. Shake the container to remove any lumps in the powder (fig. 1.6).
5. Wait 10 seconds before opening the container. Using the dosing spoon, collect the powder without compressing it and level off with the flat part of the spatula (fig. 1.7). Pour the collected alginate into a rubber mixing bowl.
6. Repeat the operation until you reach the desired amount.
7. Seal the container properly at the end of dosing (fig. 1.8).
8. Add water to the powder. Pour 1/3 cup of water (18 ml) for each spoonful of powder (9 g). Example: for 2 measures of powder (18 g) sufficient for an average upper impression tray, pour 2/3 cup of water (36 ml) (fig. 1.9).
9. With the aid of a spatula, mix the compound to achieve an even consistency (fig. 1.10), observing the times stated in the technical data table (see paragraph 14).

B. 18 g SINGLE-DOSE BAG (Orthoprint/ Tropicalgin)

1. Choose a suitable impression tray.
2. Open the single-dose bag (fig. 2.1).
3. Completely empty the content of the bag (18 g) into a rubber mixing bowl (fig. 2.2).
4. Add 36 ml of water.
5. With the aid of a spatula, mix the compound to achieve an even consistency (fig. 2.3), observing the times stated in the technical data table (see paragraph 14).

12.3 IMPRESSION TAKING PROCEDURE

1. Place a suitable amount of material on the impression tray (fig. 3.1).
2. Place the loaded impression tray inside the patient's mouth within the working time (see paragraph 14, technical data table) (fig. 3.2).
3. Remove the impression from the patient's mouth once the setting time has been reached (see paragraph 14, technical data table).
4. Check that the impression is intact, have the patient rinse their mouth and check there are no material residues to be removed.

12.4 MANUAL IMPRESSION CLEANING AND DISINFECTION

WARNINGS/PRECAUTIONS REGARDING MANUAL IMPRESSION CLEANING AND DISINFECTION

- Always disinfect impressions manually with a specific disinfectant for alginate, using either a concentrated product containing quaternary ammonium salts and synergistic co-formulants (formulation for immersion), or a ready-to-use product containing mixtures of alcohol and surface tension reducers (spray formulation). Follow the chosen disinfectant manufacturer's instructions for contact time.
- Using an unsuitable disinfectant or inappropriate solution may compromise the impression.
- Do not use a toothbrush to remove dirt as this may damage the impression.
- Manual cleaning alone is not sufficient for proper impression treatment, and the manual disinfection step must always be carried out.
- Impressions MUST NOT BE:
 - o sterilised in an autoclave, with chemical vapour, in dry heat or with cold chemical immersion sterilisers;
 - o cleaned/disinfected in instrument washing machines or thermal disinfection devices;
 - o cleaned in ultrasound baths.
- No automatic methods of cleaning and disinfecting impressions have been validated.
- Impressions must always be cleaned and disinfected before being sent to the laboratory for subsequent processing.
- The appropriately cleaned and disinfected impression must be stored in a closed, clean place protected against splashes, spray or aerosols of body fluids at

room temperatures of approximately 23°C / 73°F.

12.4.1 PREPARATION FOR CLEANING AND DISINFECTION

1. Put on a pair of disposable gloves, facemask, goggles and work clothes.
2. If the impression is to be disinfected by immersion (see paragraph 12.4.2 A), prepare a tray containing a solution of impression cleaner/disinfectant that has been approved pursuant to local regulations, containing quaternary ammonium compounds and synergistic co-formulants (e.g. Zeta 7 Solution).

12.4.2 MANUAL CLEANING AND DISINFECTION OF IMPRESSIONS

A. IMMERSION (fig. 4.2 A)

Use the immersion cleaner-disinfectant following the manufacturer's instructions for use.

1. Rinse the impression under running water **for 30 seconds** immediately after removing it from the patient's mouth (fig. 4.1). Carry out the cleaning step within 5 minutes of removal from the mouth. If necessary, increase the rinse time until all visible residues have been removed.
2. Immerse the impression in the cleaner-disinfectant solution (prepared previously as described in point 2 of paragraph 12.4.1. PREPARATION FOR CLEANING AND DISINFECTION).
3. Leave the impression in the solution for the time indicated in the instructions for use of the cleaner-disinfectant.
4. Rinse and remove excess water.

B. SPRAY (fig. 4.2 B)

Use the cleaner-disinfectant spray following the manufacturer's instructions for use.

1. Rinse the impression under running water **for 60 seconds** immediately after removing it from the patient's mouth (fig. 4.1). Carry out the cleaning step within 5 minutes of removal from the mouth. If necessary, increase the rinse time until all visible residues have been removed.
2. Place the impression inside a transparent bag.
3. Introduce the nozzle of the bottle of impression cleaner/disinfectant spray, approved in accordance with local regulations, containing ethanol, 2-propanol and co-formulants (e.g. Zeta 7 Spray), into the transparent bag.
4. Hold the bag closed to prevent the aerosol escaping.
5. Spray the cleaner-disinfectant, taking care to cover the whole surface of the impression and impression tray and observe the contact time stated in the instructions for use of the cleaner-disinfectant.
6. Hold the bag closed and leave to act.
7. Remove the impression from the bag and dispose of the bag.
8. Rinse and remove excess water.

13. INFORMATION ON THE PROPER MANAGEMENT OF PACKAGING AND ACCESSORIES

Use clean and appropriately disinfected and sterilised or ready-to-use accessories.

The packaging and accessories listed under point 13.1 can be reused safely without being cleaned and disinfected, provided the following guidelines for conservation, dosing and storage of components are observed:

- Handle, dispense and store the packaging and accessories outside the surgical area in a clean place away from splashes, sprays or aerosols of body fluids;
- Always wear new, uncontaminated gloves when handling packaging and accessories. If gloves become contaminated, change them immediately before handling the product;
- Close packaging immediately after dosing.

If contaminated, dispose of products and accessories as described in paragraph 15 on disposal.

13.1 PACKAGING AND ACCESSORIES THAT DO NOT REQUIRE CLEANING AND DISINFECTION

- 453 g/500 g bag
- Container
- Dosing set

14. TECHNICAL DATA

	Hydrogum	Orthoprint	Tropicalgin Normal	Tropicalgin	Tropicalgin Extra Fast	Neocolloid
Powder-to-water mixing ratio	9 g:18 ml	9 g:18 ml	9 g:18 ml	9 g:18 ml	9 g:18 ml	9 g:18 ml
Mixing time	30"	30"	45"	45"	30"	45"
Working time* (including mixing time)	1'10"	1'05"	1'50"	1'35"	1'15"	2'00"
Time in mouth (min.)	1'00"	45"	1'00"	1'00"	25"	1'30"
Setting time*	2'10"	1'50"	2'50"	2'35"	1'40"	3'30"
Elastic recovery	96.5%	96.5%	97%	97.5%	97%	97%
Tear resistance (N/mm)	0.7	0.8	0.7	0.7	0.7	0.65
Compression strain	11.5%	13%	13%	11%	13%	12.5%
Initial setting time (23°C/73°F)	2'20"	1'55"	2'30"	2'15"	1'30"	3'30"

The water temperature and hardness may affect the times for clinical use of the material. Cold water retards curing; hot water accelerates it.

*The times mentioned above are intended from the start of mixing at 23°C/73°F. Higher temperatures reduce these times, lower temperatures extend them.

15. STORAGE, STABILITY AND DISPOSAL

The devices should be used within dental surgeries (public or private health facilities with the necessary legal authorisations).

Based on the information provided in the SDS, the waste generated by disposal of the product can cause physical harm to people and the environment. It should therefore be managed as hazardous waste according to local regulations. If contaminated, the devices and the impression must be disposed of as special waste at risk of biological contamination. For instructions regarding handling of contaminated components see paragraph 12.4. Always wear gloves when handling the device. If the devices and impression are not contaminated, dispose of them in accordance with applicable regulations.

The expiry date relates to the unopened package properly stored. Store the product in an airtight container. Close the container at the end of dosing and keep it at a temperature of between 5 and 27°C / 41-80°F, in a dry place protected from direct sunlight.

16. IMPRESSION CASTING

- Hydrogum, Orthoprint and Neocolloid: it is recommended to cast the gypsum immediately after disinfection or within 48 hours of impression-taking.
- Tropicalgin, Tropicalgin Normal: when products containing quaternary ammonium compounds and synergistic co-formulants are used for disinfection, it is recommended to cast the gypsum immediately after disinfection or within 24 hours of impression-taking. When products containing mixtures of alcohol and surface tension reducers are used for disinfection, cast the gypsum within 5 hours.
- Tropicalgin Extra Fast: when products containing quaternary ammonium compounds and synergistic co-formulants are used for disinfection, it is recommended to cast the gypsum immediately after disinfection or within 20 hours of impression-taking. When products containing mixtures of alcohol and surface tension reducers are used for disinfection, cast the gypsum within 5 hours.
- The material is compatible with study model or working model gypsums (type 3 gypsum e.g. Elite Model - Zhermack or type 4 gypsums e.g. Elite Rock - Zhermack).

17. CONSERVATION OF THE IMPRESSION

If the gypsum cannot be cast immediately, remove the water in excess and store the impression at room temperature (23°C) in a clean polyethylene bag with an airtight seal.

18. IMPORTANT REMARKS

Information provided in any way, even during demos, does not invalidate the instructions for use. Operators are required to check that the product is suitable for the application envisaged. The manufacturer cannot be held responsible for damage, including to third parties, deriving from failure to follow instructions or from unsuitability for an application. The manufacturer's liability is, in any case, limited to the value of the products supplied. Report any serious incident involving the medical device to the manufacturer and to the relevant authorities.

453 g/500 g BAG (Hydrogum/ Orthoprint/ Tropicalgin/ Tropicalgin Normal/ Tropicalgin Extra Fast/ Neocolloid)

Five-year shelf-life from the production date.

18g SINGLE-DOSE BAG (Orthoprint/ Tropicalgin)

Three-year shelf-life from the production date.

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