# Occlufast Rock

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
Occlufast Rock 90±5 Shore A


[dentistrysirona.com/zhermack](http://dentistrysirona.com/zhermack)

**PREPARING THE DEVICES: STEP-BY-STEP INSTRUCTIONS**

**DISINFECTING AND CLEANING THE DEVICE: GENERAL STEP-BY-STEP INSTRUCTIONS**

**ATTENTION:** For correct manual cleaning and disinfection of the cartridge and mixing tip read sections 9.5.2 and 9.5.3 carefully

## EN – Occlufast Rock

### Instructions for use

#### 1. INTENDED USE

A-silicone bite registration material.

#### 2. PRODUCT DESCRIPTION

Occlufast Rock has a hardness of  $90 \pm 5$  Sh A.

#### 3. PACKAGING

50 ml cartridge.

#### 4. COMPOSITION

Vinyl polysiloxane, silica fillers, inorganic aluminium fillers, methylhydrogenosiloxane, organoplatinum complex, pigments.

#### 5. INDICATIONS FOR USE

The devices are intended to be used by dental professionals (user licensed to operate on the patient) for taking bite registration on patients from 3 years of age in a professional dental environment (public or private health facilities with the necessary regulatory approvals).

#### 6. CONTRAINDICATIONS

Do not use on patients with known hypersensitivity to any of the components. For more information on any allergens contained in the scents, where present, please contact Zhermack.

#### 7. SIDE EFFECTS

Irritation, redness or signs of hypersensitivity may occur in case of allergy to any of the components.

#### 8. CLINICAL BENEFITS

Dimensionally stable material permits intraoral registration of the relative positions of the arches.

#### GENERAL PRECAUTIONS/WARNINGS:

- The product's instructions for use must be kept for the duration of its use.
- The devices should be used by professionals operating within a professional dental environment (legally-certified public or private health facilities).
- The hardening of Addition Silicones (Vinyl polysiloxanes) is inhibited by contact with latex. **Avoid direct or indirect contact between Addition Silicones (Vinyl polysiloxanes) and gloves or other items made of latex. In the event of accidental contact with latex**, wash hands thoroughly to eliminate all traces of contaminants and use vinyl or nitrile gloves.
- In order to reduce the risk of cross-contamination, always use new, uncontaminated gloves when handling cartridges and accessories.
- Do not use the material:
  - if the batch number and expiry date on the outer packaging are missing or illegible;
  - if the primary packaging (cartridge) is damaged.
- Safety Data Sheet available subject to request.

#### 9. STEP-BY-STEP INSTRUCTIONS

##### 9.1 PRELIMINARY OPERATIONS

Wear appropriate protective goggles, face mask, work clothes and gloves. Protective goggles are also recommended for the patients.

##### 9.2 DEVICE ASSEMBLY

#### PRECAUTIONS/ WARNINGS CONCERNING DEVICE ASSEMBLY:

- Do not use after the expiry date indicated on the packaging.
  - The product has not been tested for compatibility with accessories other than those with which it is supplied. Use of non-original accessories may cause one or more sequences of events that could lead to a reduction in the precision of bite registration. The use of non-original accessories is therefore not recommended.
  - In order to reduce the risk of cross-contamination, before use on each new patient, always apply a disposable cover/barrier system to the dispenser and cartridge, to be removed after each use. The disposable cover/ barrier system does not replace the recommended cleaning and disinfection operations.
1. Push the black lever underneath the rack upwards and pull the dispenser rack as far as it will go (fig. 1.1).
  2. Lift the cartridge lock (fig. 1.2).
  3. Insert the cartridge making sure that it fits the right direction. The cartridge must be inserted with the "V" on its base facing downwards (fig. 1.2). The cartridge lock can now be lowered easily.
  4. After placing the cartridge into the dispenser, turn the cap by a 1/4 turn, pull it out with a firm movement and throw it away (fig. 1.3).
  5. Press the dispenser lever gently to squeeze out a small amount of the two components until they both emerge at the same time (fig. 1.4). **Clean by removing all excess product, taking care not to contaminate the two outlet holes.**
  6. Select the most appropriate mixing tip; the colour of the tip must be the same as the colour of the cartridge cap.
  7. Insert the mixing tip in a clean motion so that the V-shaped notch on the coloured ring aligns with its housing on the cartridge (fig. 1.5). Turn the tip clockwise by a 1/4 turn, exerting pressure only on the coloured ring and not on the mixing tip itself (fig. 1.6).
  8. Exert a slight pressure on the dispenser lever to extrude a sufficient quantity of material. Do not use if the mixture is not smooth. In this case, remove the dispenser from the operating area and eliminate the mixing tip, extrude a small quantity of material to level off the component and insert a new mixing tip.
  9. Before using on a new patient, always apply a disposable cover/barrier system to the cartridge, to be removed after every use (fig. 1.7). The cover/barrier system is not a substitute for the recommended cleaning and disinfection operations.
  10. Proceed to use the device on the patient (see paragraph 9.3).

### 9.3 USING THE DEVICE ON THE PATIENT

#### **WARNING/PRECAUTIONS REGARDING USE OF THE DEVICE ON THE PATIENT:**

- The product is intended for contact with intact mucosae.
- Use the product at an ambient temperature of 23°C/73°F (higher temperatures reduce working time, lower temperatures extend working time).
- Always (in any case before each use) level the material (fig. 1.4) before attaching the mixing tip. When applying the material, take care to keep the mixing tip immersed inside it.
- Use on patients sensitive to Addition Silicones (Vinyl polysiloxanes) is not recommended. If irritation, redness or other signs of hypersensitivity occur, stop using the product immediately and take the steps necessary to ensure the patient's safety.
- In order to avoid complications, it is advisable to secure undercuts before bite registration.
- The hardening of Addition Silicones (Vinyl polysiloxanes) is inhibited by contact with latex, sulphide-based gingival retraction solutions, polyethers and the unpolymerised surface of acrylic or bis-acryl resins. Avoid direct or indirect contact of the material with these substances. Should the patient's hard and soft oral tissues or accessory surfaces with which the material comes into contact be contaminated with these substances, clean them thoroughly with water, or with ethyl alcohol/other solvents in the case of acrylic or bis-acryl resins.
- Avoid contact between the cartridge and the patient's mouth.

1. Secure undercuts before using the product.
2. Press the dispenser lever gently to squeeze out an appropriate amount of material. Do not use the material if it is visibly inconsistent.
3. Apply the product using the desired technique (occlusal or vestibular).
4. When using the occlusal technique, having applied the material to the occlusal surface of the teeth, the dentist will ask the patient to bite within the working time indicated in table 10 Technical Data.
5. When using the vestibular technique, the product must be applied to the teeth from the vestibular side, with the dental arches in the position to be registered.
6. The material will not be removed from the patient's mouth until the end of the time in mouth indicated in table 10 Technical Data.
7. Check that the bite registration is intact and make sure there are no residues to be removed from the patient's mouth.
8. Proceed with the cleaning and disinfection phases described in **sections 9.5.2 and 9.5.3.**

### 9.4 MANUAL CLEANING AND DISINFECTION OF THE BITE REGISTRATION

#### **WARNINGS/ PRECAUTIONS REGARDING THE MANUAL CLEANING AND DISINFECTION OF THE BITE REGISTRATION:**

- Manual cleaning alone is not sufficient for correct treatment of the bite registration and the manual disinfection phase should always be performed as well.
- Always disinfect bite registrations manually using a specific disinfectant for addition silicones (vinyl polysiloxanes) for immersion (concentrated disinfectant containing quaternary ammonium salts and synergistic co-formulant) or by spraying a ready-to-use disinfectant containing alcohol mixtures and surface tension reducers. Follow the chosen disinfectant manufacturer's instructions for the contact times.
- Store the clean and disinfected bite registration in a closed, dry and clean place, away from splashes, jets or aerosols of bodily fluids, at an ambient temperature of around 23°C / 73°F.

#### 9.4.1 PREPARATION BEFORE THE CLEANING AND DISINFECTION PHASES

1. Wear a pair of disposable gloves, face mask, goggles and work clothes.
2. For disinfection of the bite registration through immersion (see paragraph **9.4.3 A**), prepare a specific cleaning and disinfecting solution for bites, approved according to local regulations, containing quaternary ammonium salt compounds and synergistic coformulants (e.g. Zeta 7 Solution).

#### 9.4.2 MANUAL CLEANING

1. Rinse the bite registration under running water **for 30 seconds** within a maximum of 5 minutes after removal from the mouth (fig. 2.1). If necessary, increase the rinsing time until no residues are visible.
2. Proceed with disinfection through immersion (see paragraph **9.4.3 A**) or spraying (see paragraph **9.4.3 B**).

#### 9.4.3 MANUAL DISINFECTION OF BITE REGISTRATIONS

##### **A. IMMERSION (fig. 2.2 A)**

1. After the cleaning procedure (see paragraph 9.4.2), immerse the bite registration in the cleaning and disinfecting solution (see point 2 of paragraph 9.4.1.) for the envisaged contact time.
2. Proceed with rinsing and/or drying where envisaged.

##### **B. SPRAY (fig. 2.2 B)**

1. After the cleaning procedure (see paragraph 9.4.2), place the bite registration in a transparent bag.
2. Introduce into the bag the nozzle of the bottle of specific cleaner-disinfectant spray for bites, approved according to local regulations, containing ethanol, 2-propanol and coformulants (e.g. Zeta 7 Spray).
3. Keeping the bag closed, spray the cleaner-disinfectant onto all the surfaces of the bite registration, complying with the cleaner-disinfectant contact times.
4. Keep the bag closed and leave to act.
5. Extract the bite registration from the bag and dispose of the bag.
6. Proceed with rinsing and/or drying where envisaged.

### 9.5 MANUAL CLEANING AND DISINFECTION OF THE CARTRIDGE AND MIXING TIP

#### **WARNINGS/ PRECAUTIONS CONCERNING MANUAL CLEANING AND DISINFECTION OF THE CARTRIDGE AND MIXING TIP:**

- The protective barriers are designed for single use and must be disposed of after use, in accordance with local regulations.
- The cartridges and mixing tips must be cleaned and disinfected as described below before storage. Manual cleaning alone is not sufficient for correct treatment of the cartridge with the mixing tip inserted. Keep the used, clean and disinfected mixing tip on the cartridge as a cap.

- Use only a disinfectant solution that is approved for its efficacy and use in accordance with the IFU of the disinfectant solution manufacturer.
- Perform manual cleaning and disinfection of the cartridge with the mixing tip attached within 1 hour of use.
- Remove the used mixing tip from the cartridge before reuse and discard it. Level off the cartridge and insert a new mixing tip.
- The end of the cartridge life cycle is normally determined by the use-by date and/or consumption of the material it contains.

#### 9.5.1 DISASSEMBLY OF THE DEVICE AFTER USE

1. In the area used for cleaning and disinfection, remove the intra-oral tip of the mixing tip and any residues of bite material present inside it.
2. Remove the disposable cover/ barrier system.
3. Remove the contaminated gloves, sanitise hands and put on a clean pair of gloves.
4. Remove the cartridge with the mixing tip attached from the dispenser, taking care to avoid direct contact with the contaminated tip.
5. Clean and disinfect the dispenser according to the dispenser manufacturer's instructions.

#### 9.5.2 MANUAL CLEANING

1. Do not remove the mixing tip from the cartridge.
2. Clean using disposable wipes completely impregnated with a cleaning and disinfecting solution\*.
  - a. Take a disposable wipe completely impregnated with a cleaning and disinfecting solution\* and wipe it **twice over all surfaces of the cartridge**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.
  - b. Take a new disposable wipe completely impregnated with a cleaning and disinfecting solution\* and wipe it **twice over all surfaces of the mixing tip of the cartridge**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.

Use additional wipes if necessary.

3. Use a clean soft brush (e.g. Interlock® REF 09084) wetted with a cleaning and disinfecting solution\* to clean the parts of the cartridge and the mixing tip with recesses and/or grooves (e.g.: connection between the mixing tip and the head of the cartridge) taking care to remove all visible residues. Do not scrub the cartridge label with the brush.
4. Remove the residues of cleaning and disinfecting solution\* with a wet cloth. Proceed with disinfection.

#### 9.5.3 MANUAL DISINFECTION

1. Do not remove the mixing tip from the cartridge.
2. Disinfect using disposable wipes completely impregnated with a cleaning and disinfecting solution\*.
  - a. Take a disposable wipe completely impregnated with a cleaning and disinfecting solution\* and wipe it **twice over all surfaces of the cartridge**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.
  - b. Take a new disposable wipe completely impregnated with a cleaning and disinfecting solution\* and wipe it **twice over all surfaces of the mixing tip of the cartridge**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.

Use additional wipes if necessary.

3. Allow the device to dry naturally for 5 minutes.
4. Inspect the cartridge and tip to make sure that all visible contamination has been removed. If visible contamination is still present, repeat the cleaning and disinfection steps described in points 9.5.2 and 9.5.3.
5. Before reusing the cartridge, insert it into the dispenser again, remove and discard the used mixing tip (see paragraph 9.2).

The instructions provided above have been validated by the manufacturer of the medical device as the provider of a reusable medical device. However, the responsibility for correct performance of the cleaning and disinfection phases of the medical device lies with the person who performs the process, who must guarantee that it is carried out correctly in the facility where the cleaning and disinfection phases take place and that the desired result is achieved. This requires verification and/or validation and continuous monitoring of the process.

\*cleaning and disinfecting solution, bactericide, yeasticide, tuberculocide and virucide containing ethyl alcohol (~35%) and isopropyl alcohol (~15%) mixtures (e.g. Zeta 3 Soft, cleaning and disinfecting solution, 5-minute contact time) approved in accordance with local regulations; use in accordance with the instructions for use of the disinfectant solution manufacturer.

### 10. TECHNICAL DATA

|                                                 | Occlufast Rock                     |
|-------------------------------------------------|------------------------------------|
| ISO 4823                                        | Type B Bite registration materials |
| Working time*                                   | 30"                                |
| Time in mouth**                                 | 1'15"                              |
| Shore A hardness 1 hour                         | 90±5                               |
| Detail reproduction                             | 20 µm                              |
| Linear dimensional variation (approx.) 24 hours | 0,25%                              |

\* The time is intended from the start of mixing at 23°C/73°F. Higher temperatures reduce this time, lower temperatures extend it.

\*\* The time in mouth is intended as at 35°C / 95°F.

### 11. STORAGE, STABILITY AND DISPOSAL

The devices are intended for use within a professional dental environment (public-sector or private-sector healthcare facilities possessing the necessary regulatory permits).

- Dispose of the cartridge when empty. Do not attempt to reload used cartridges. The intra-oral tip with the material inside it must be disposed of as waste with a risk of biological contamination.

• Based on the information contained in the SDS, the waste produced by disposal of the product does not present any physical hazards for people or for the environment. It may therefore be managed in accordance with applicable local regulations, as non-hazardous waste. For handling of contaminated components, see paragraphs 9.4 and 9.5. If the material and bite are contaminated, dispose of them as special waste with a risk of biological contamination. If they are not contaminated, dispose of them in accordance with applicable local regulations.

Do not use after the expiry date indicated on the packaging. Store the product at a temperature between 5°C/41°F and 27°C/ 80°F, away from direct sunlight.

## **12. STORAGE OF THE BITE REGISTRATION**

When disinfection is complete, dry the registration and use it within a maximum of 7 days. The registration can be sent to the laboratory as a normal bite, without particular precautions, provided it is protected against mechanical stress.

## **13. IMPORTANT REMARKS**

Information provided in any way, even during demonstrations, does not invalidate the instructions for use. Operators are required to check that the product is suitable for the application envisaged. The manufacturer cannot be held responsible for damage, including to third parties, deriving from failure to follow instructions or from unsuitability for an application. The manufacturer's liability is, in any case, limited to the value of the products supplied. Report any serious incident involving the medical device to the manufacturer and to the relevant authorities.

Three-year shelf-life from the production date.

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